



Email : sales@laboratoryinstrumentindia.com

Product Name :
1.6l Engine Bench

Product Code :
LBNY-0005-121000284



Description :

Trainer is mounted on a mobile stand constructed of 2(5cm) tubular steel equipped with four 5(12.5cm) phenol casters (2-fixed and 2-swivel w/locks). 1.6l Engine Bench Trainer is equipped with safety guards covering all rotating components. The engine performance trainer is a fully operational engine mounted on a mobile stand. The trainer is easy to work on and all major components, wiring and connectors match OEM configuration. Trainer is based on 2007 and newer 1.6L global engine platform.

Advantages:

Engine systems respond to inserted faults with real world symptoms, OEM DTCs and check engine light operation.

Demonstration of systems to groups of students without the access limitations of a complete vehicle.

Features:

Sequential fuel injection w/coil-on-plug

1.6L DOHC engine

Vehicle OEM VIN number supplied with trainer

Adaptive charging system

Technical Specification :

Fuel requirements: 87 octane gasoline

Dimensions: 29 W x 48 L x 50 H (72.5 x 120 x 125 cm) w/o pallet 35 W x 60 L x 56 H (90 x 150 x 140 cm) w/pallet

Weight: 775 lbs (340.9 kg) w/o pallet 850 lbs (386.4 kg) w/pallet.



Laboratory Instrument India