



Email : sales@laboratoryinstrumentindia.com

Product Name :
2.2L Swivel Engine Stand Starting Station

Product Code :
LBNY-0005-121000288



Description :

The 2.2L Swivel Engine Stand Starting Station for student instruction and practice covering engine mechanical inspection, disassembly, measurement and reassembly. A steel mounting plate is used for engine mounting that allows full access for rear main seal and flywheel service. Each engine comes equipped with a manual transmission flywheel. The engine is securely mounted to a heavy duty rotating device equipped with an automatic lock to allow full 360 degree rotation. All necessary engine systems connections are provided to allow live operation of each swivel engine stand. The engine starting station contains all components and systems necessary for the starting and running of the 2.2L engine.

Advantages:

Provides a single starting station that can be used with multiple engine stands
All rotating components are covered with safety shields or guards to prevent injuries
Accessibility of engine components is enhanced by not having the requirement of whole vehicles in the shop
Ability to verify proper engine systems operation before and after engine repair
Continuity of subject matter is provided by students working on identical engines.

Technical Specification :

Dimensions: 56 W x 40 L x 56 H w/o pallet / 59 W x 46 L x 59 H w/pallet
Weight: 800 lbs w/o pallet / 900 lbs w/pallet.

Laboratory Instrument India

