



Email : sales@laboratoryinstrumentindia.com

Product Name :

Workshop Universal Tensile Material Testing Machine

Product Code :

LBNY-0004-730008



Description :

Workshop Universal Tensile Material Testing Machine

Technical Specification :

Universal Tensile Material Testing Machine

Power = 1.5KW

Testing Accuracy = 0.5

Size = 760x530x1830(mm)

Machine Weight = 330KG

Test value = 50KN

Effective test width = 400mm

Displacement resolution = 0.01mm

Displacement speed adjustment range = 0.01-500mm/min

Displacement indication limit error = Within $\pm 1\%$ of indication

Test force measurement range = 2% to 100% F.S

Workshop Universal Tensile Material Testing Machine, Workshop Universal Tensile Material Testing Machine Equipments, Workshop Universal Tensile Material Testing Machine Tools, Workshop Universal Tensile Material Testing Machine Tool Kits, Workshop Universal Tensile Material Testing Machine Manufacturers, Workshop Universal Tensile Material Testing Machine Suppliers from India, China, Kenya



Laboratory Instrument India