



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cutaway Turbocharger

**Product Code :**  
LBNY-0005-121000290



**Description :**

The Cutaway Turbocharger is used to demonstrate the operation of a variable displacement turbocharger. The Cutaway Turbocharger trainer is cutaway in such a manner as to clearly show the operation of internal parts that provide the variable displacement. This component utilizes a variable output displacement by using variable vane technology (VVT).

**Advantages:**

All internal components are visible for better understanding of turbo operation  
Students can observe the movement of the variable-vane turbine  
Fully functional vacuum diaphragm that moves turbo plate.

**Main Features:**

Operational vacuum diaphragm  
Variable vane technology (VVT)  
Clear demonstration of all internal turbocharger components.

**Technical Specification :**

Cutaway Turbocharger

---

# Laboratory Instrument India

