



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Brinell Hardness Testing Device

Product Code :
LBNY-0004-730009



Description :

Workshop Brinell Hardness Testing Device

Technical Specification :

Brinell Hardness Testing Device

Automatic testing after preliminary test force is loaded by hand. Measure indentation by measuring microscope then find the hardness value from table. Electronic close-loop load control system.

Test force: 62.5kgf(612.9N)100kgf(980.7N)125kgf(1226N)187.5kgf(1839N)

Measurement unit: 0.005mm

Hardness test range: 8~650HBW

Executive standard: GB/T231.2,JJG150

Indenter - Upper wall distance: 55mm

Brinell ruler:HBW2.5/187.5 HBW5/125 HBW5/750 HBW10/100 HBW10/3

Height of test piece: 220mm

Burden time:0~60s

Total magnification: 20X

Indenter - outside wall distance:135mm

Workshop Brinell Hardness Testing Device, Workshop Brinell Hardness Testing Device Equipments, Workshop Brinell Hardness Testing Device Tools, Workshop Brinell Hardness Testing Device Tool Kits, Workshop Brinell Hardness Testing Device Manufacturers, Workshop Brinell Hardness Testing Device Suppliers from India, China, Kenya



Laboratory Instrument India