



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Electric Power Steering Trainer

**Product Code :**  
LBNY-0005-121000302



#### **Description :**

The Electric Power Steering Trainer provides a safe and effective method of demonstrating the principles of electric power steering systems used on today's vehicles. The Electric Power Steering Trainer is a fully functional electric power steering trainer.

#### **Advantages:**

Provides easy access to all system components without the confinements of a real vehicle  
Demonstrates full functionality of a typical electric power steering system  
Instructor can manually insert up to four system faults for student troubleshooting and diagnosis training and practice.  
Students can feel the power steering assist effort being applied.

#### **Features:**

Fully functional system with live operation  
Built in ammeter showing assist torque applied  
Requires a 12V automotive battery as a power source  
All OEM parts are used thus making instruction relative to today's vehicles

#### **Technical Specification :**

Dimensions: 22 L x 22 W x 9 H (55.9 x 55.9 x 22.8 cm)  
Weight: 15 lb (6.8 kg).



**Laboratory Instrument India**