



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electric Power Steering Trainer

Product Code :
LBNY-0005-121000302



Description :

The Electric Power Steering Trainer provides a safe and effective method of demonstrating the principles of electric power steering systems used on today's vehicles. The Electric Power Steering Trainer is a fully functional electric power steering trainer.

Advantages:

Provides easy access to all system components without the confinements of a real vehicle
Demonstrates full functionality of a typical electric power steering system
Instructor can manually insert up to four system faults for student troubleshooting and diagnosis training and practice.
Students can feel the power steering assist effort being applied.

Features:

Fully functional system with live operation
Built in ammeter showing assist torque applied
Requires a 12V automotive battery as a power source
All OEM parts are used thus making instruction relative to today's vehicles

Technical Specification :

Dimensions: 22 L x 22 W x 9 H (55.9 x 55.9 x 22.8 cm)
Weight: 15 lb (6.8 kg).



Laboratory Instrument India