



Email : sales@laboratoryinstrumentindia.com

Product Name :
Overdrive Cutaway

Product Code :
LBNY-0005-121000305



Description :

Epicyclic overdrive cutaway model, with electromagnetic actuator, of the type fitted to Volvo motorcars.

Technical Specification :

Equipment Layout: On Base
Operated: Manually (Through a Crank Handle)
Dimensions: 50 x 30 x 40(h) cm
Net Weight: 15 kg
Gross Weight: 20 kg.



Laboratory Instrument India