



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Charpy Pendulum Impact Testing Machine

**Product Code :**  
LBNY-0004-730010



**Description :**

Workshop Charpy Pendulum Impact Testing Machine

**Technical Specification :**

Charpy Pendulum Impact Testing Machine

Power supply: 220V 50Hz

Application: charpy impact tester

Name: impact testing machines

Function: pendulum impact testing machine

Impact speed = 2.9m / s or 3.8M / S

Distance from pendulum center to sample center = 221mm, 380mm

Pendulum angle = 150°

Jaw support spacing = 40mm, 60mm, 70mm, 95mm

Support blade corner radius = R = 1mm

Impact blade angle = 30 °

Impact blade fillet radius = R = 2mm ± 0.5mm:

Workshop Charpy Pendulum Impact Testing Machine, Workshop Charpy Pendulum Impact Testing Machine Equipments, Workshop Charpy Pendulum Impact Testing Machine Tools, Workshop Charpy Pendulum Impact Testing Machine Tool Kits, Workshop Charpy Pendulum Impact Testing Machine Manufacturers, Workshop Charpy Pendulum Impact Testing Machine Suppliers from India, China, Kenya



**Laboratory Instrument India**