



Email : sales@laboratoryinstrumentindia.com

Product Name :
Self Locking Differential Cutaway

Product Code :
LBNY-0005-121000311



www.laboratoryinstrumentindia.com

Description :

Accurate section of an assembly including hypoid bevel gear pair with self-locking differential.

Technical Specification :

Equipment Layout: On Table Stand
Operated: Manually (Through a Crank Handle)
Dimensions: 40 x 30 x 40(h) cm
Net Weight: 25 kg
Gross Weight: 30 kg.



Laboratory Instrument India