



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Self Locking Differential Cutaway

**Product Code :**  
LBNY-0005-121000311



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Accurate section of an assembly including hypoid bevel gear pair with self-locking differential.

**Technical Specification :**

Equipment Layout: On Table Stand  
Operated: Manually (Through a Crank Handle)  
Dimensions: 40 x 30 x 40(h) cm  
Net Weight: 25 kg  
Gross Weight: 30 kg.



**Laboratory Instrument India**