



Email : sales@laboratoryinstrumentindia.com

Product Name :
Four WD Vehicle Transmission Assembly With Five Speed
Manual Gearbox Cutaway

Product Code :
LBNY-0005-121000316



Description :

The various parts are supplied as accurate real cutaway, connected together to show their operation in an easy and immediately understandable manner. Transmission assembly fitted to most 4WD vehicles.

Technical Specification :

2-Speeds reduction gear
Movement restorer with front wheel drive manual control
Gearbox: 5 speeds and reverse
Self-locking hypoid differentials
Drive shafts with universal joints
Manual operation.

Other Specifications:

Equipment Layout: On Stand with Wheels
Operated: Manually (Through a Crank Handle)
Dimensions: 75 x 225 x 105(h) cm
Net Weight: 230 kg
Gross Weight: 300 kg.

Laboratory Instrument India

