



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Four WD Transfer Case With Locking Differential Cutaway

**Product Code :**  
LBNY-0005-121000325



**Description :**

It is permanent transfer with possibility of differential locking in order to make the front axle independent from the rear one. Speed reducer transfer case for 4WD cars.

**Technical Specification :**

Equipment Layout: On Base  
Operated: Manually (Through a Crank Handle)  
Dimensions: 40 x 40 x 60(h) cm  
Net Weight: 25 kg  
Gross Weight: 40 kg.



**Laboratory Instrument India**