



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop MMA (Manual Metal Arc) Machine

Product Code :
LBNY-0004-740001



Description :

Workshop MMA (Manual Metal Arc) Machine, Workshop MMA (Manual Metal Arc) Machine Equipments, Workshop MMA (Manual Metal Arc) Machine Tools, Workshop MMA (Manual Metal Arc) Machine Tool Kits, Workshop MMA (Manual Metal Arc) Machine Manufacturers, Workshop MMA (Manual Metal Arc) Machine Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

MMA (Manual Metal Arc) Machine

Output Current DC 250 -300A

Lightweight, Suitable for decoration and maintenance industries.

It has strong anti power grid fluctuation ability and wide adaptability.

Constant power control and stable output.

It is suitable for all kinds of electrode, undervoltage, overvoltage, overcurrent and automatic protection.

IGBT inverter technology, low power consumption, low noise, high frequency, high efficiency, good antiinterference performance and stable arcing.



Laboratory Instrument India