




Email : sales@laboratoryinstrumentindia.com

Product Name :
HGV Thirteen Speed Manual Fuller Cutaway

Product Code :
LBNY-0005-121000331



Description :

HGV Thirteen Speed Manual Fuller Cutaway is a non-synchronized gearbox and it is divided in 2 parts: In the other part of the gearbox (on the output shaft side), there are the standard, low ratio and semi-low ratio speed gears, pneumatically controlled by the pre-selectors provided on the gear lever. On the engine side are the 1st 2nd 3rd and 4th speed gears, reverse speed gears and extra low ratio pick-up speed gears. All these gears are mechanically controlled by the change gear lever. Quick change gear box used in 300/400HP heavyweight vehicles with mechanical and pneumatic control.

Technical Specification :

Equipment Layout: On Stand with Wheels
Operated: Manually (Through a Crank Handle)
Dimensions: 110 x 100 x 130(h) cm
Net Weight: 300 kg
Gross Weight: 400 kg.



Laboratory Instrument India

