



Email : sales@laboratoryinstrumentindia.com

Product Name :
HGV Manual Gearbox Cutaway

Product Code :
LBNY-0005-121000335



Description :

Exposed components include: 8 speed and crawler, hypercyclic box, reverse gear, synchro hubs, selector rails, gear shafts and oil system. Fully sectioned HGV gearbox, showing the internal workings of a sectioned gearbox.

Technical Specification :

Rotated via the input shaft to demonstrate the internal rotation of a speed gearbox.
Air shift system for high and low range gears.
Powder coated paint finish.
Mounted on a self contained frame and base plate, complete with castors.

Other Specifications:

Equipment Layout: On Stand with Wheels
Operated: Manually (Through a Crank Handle)
Dimensions: 100 x 120 x 130(h) cm
Net Weight: 350 kg.



Laboratory Instrument India

