



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
HGV Manual Gearbox Cutaway

**Product Code :**  
LBNY-0005-121000335



**Description :**

Exposed components include: 8 speed and crawler, hypercyclic box, reverse gear, synchro hubs, selector rails, gear shafts and oil system. Fully sectioned HGV gearbox, showing the internal workings of a sectioned gearbox.

**Technical Specification :**

Rotated via the input shaft to demonstrate the internal rotation of a speed gearbox.  
Air shift system for high and low range gears.  
Powder coated paint finish.  
Mounted on a self contained frame and base plate, complete with castors.

**Other Specifications:**

Equipment Layout: On Stand with Wheels  
Operated: Manually (Through a Crank Handle)  
Dimensions: 100 x 120 x 130(h) cm  
Net Weight: 350 kg.



**Laboratory Instrument India**

