




Email : sales@laboratoryinstrumentindia.com

Product Name :
Sequential Twin Turbocharger Cutaway

Product Code :
LBNY-0005-121000339



Description :

Accurate section of a sequential turbo showing one turbocharger for lower engine speeds and a second turbocharger at higher engine speeds.

Technical Specification :

Equipment Layout: On Base
Operated: Manually
Dimensions: 55 x 35 x 40(h) cm
Net Weight: 8 kg
Gross Weight: 15 kg.



Laboratory Instrument India