



Email : sales@laboratoryinstrumentindia.com

Product Name :
Diesel Pumping Element Cutaway

Product Code :
LBNY-0005-121000341



Description :

Diesel Pumping Element Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts. Giant model of the pumping element of a Diesel engine. It has been carefully cross sectioned to show the cylinder inlet hole and the helical spline on the piston. Many parts have been chromium-plated and galvanized for a longer life.

Technical Specification :

Equipment Layout: Standalone
Operated: Manually
Dimensions: 30 x 6 x 6(h) cm
Net Weight: 1 kg
Gross Weight: 2 kg.



Laboratory Instrument India