



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Diesel Pumping Element Cutaway

**Product Code :**  
LBNY-0005-121000341



**Description :**

Diesel Pumping Element Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts. Giant model of the pumping element of a Diesel engine. It has been carefully cross sectioned to show the cylinder inlet hole and the helical spline on the piston. Many parts have been chromium-plated and galvanized for a longer life.

**Technical Specification :**

Equipment Layout: Standalone  
Operated: Manually  
Dimensions: 30 x 6 x 6(h) cm  
Net Weight: 1 kg  
Gross Weight: 2 kg.



**Laboratory Instrument India**