




Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Pump With Radiator Cutaway

Product Code :
LBNY-0005-121000343



Description :

Many parts have been chromium-plated and galvanized for a longer life. Water Pump With Radiator Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

Technical Specification :

Equipment Layout: On Base
Operated: Manually
Dimensions: 30 x 20 x 30(h) cm
Net Weight: 3 kg
Gross Weight: 6 kg.



Laboratory Instrument India