




Email : sales@laboratoryinstrumentindia.com

Product Name :
Four In Line Cylinders Injection Pump With Pneumatic
Speed Governor Cutaway

Product Code :
LBNY-0005-121000348



Description :

It is provided with a feeding pump. Accurate section of a pump suitable for medium displacement engine (FIAT, Mercedes) with pneumatic speed governor (rock rod or acceleration rod controlled by a diaphragm connected to the suction collector).

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 20 x 30 x 30(h) cm
Net Weight: 7 kg
Gross Weight: 10 kg
Equipment Layout: On Base.



Laboratory Instrument India