



Email : sales@laboratoryinstrumentindia.com

Product Name :
Injector Pump Cutaway

Product Code :
LBNY-0005-121000349



www.laboratoryinstrumentindia.com

Description :

Accurate section of a unit injector system for commercial vehicles, where it is possible to observe:
Pumping Element
Electromagnetic Valve HD
Duster, etc.

Technical Specification :

Equipment Layout: On Base
Operated: Static
Dimensions: 40 x 20 x 10(h) cm
Net Weight: 3 kg
Gross Weight: 4 kg.



Laboratory Instrument India