



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hybrid Transmission MG Motor Generator Toyota Prius
Cutaway

Product Code :
LBNY-0005-121000351



Description :

Motor Generator 1 (MG1) effectively controls the continuously variable transmission function of the transaxle and operates as the engine starter. The Motor Generator 1 (MG1) operates as the control element for the power splitting planetary gear set. It recharges the HV battery and also supplies electrical power to drive Motor Generator 2 (MG2).

Technical Specification :

Equipment Layout: On Stand with Wheels
Operated: Manually (Through a Crank Handle)
Dimensions: 60 x 70 x 65(h) cm
Net Weight: 100 kg
Gross Weight: 130 kg.



Laboratory Instrument India