




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Direct Shift Gearbox DSG Cutaway

**Product Code :**  
LBNY-0005-121000354



**Description :**

The shafts connected to the clutches can transmit motion to an auxiliary shaft having the relevant speed-gears by means of a mechanism. Direct Shift Gearbox DSG Cutaway Also referred to as double-clutch gearbox, i Direct Shift Gearbox DSG Cutaway has two main shafts; one for the odd gears and the other for the even gears. Each shaft is inserted into one of the discs of the coaxial clutches.

**Technical Specification :**

Operated: Manually (Through a Crank Handle)  
Dimensions: 80 x 40 x 70(h) cm  
Net Weight: 60 kg  
Gross Weight: 85 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**