




Email : sales@laboratoryinstrumentindia.com

Product Name :
Vane Compressor Cutaway

Product Code :
LBNY-0005-121000357



Description :

Many parts have been chromium-plated and galvanized for a longer life. Vane Compressor Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts.

Technical Specification :

Operated: Manually
Dimensions: 20 x 20 x 20(h) cm
Net Weight: 7 kg
Gross Weight: 12 kg
Equipment Layout: On Base.



Laboratory Instrument India