



Email : sales@laboratoryinstrumentindia.com

Product Name :
High Pressure Pumps Cutaways

Product Code :
LBNY-0005-121000361



Description :

Many parts have been chromium-plated and galvanized for a longer life. High Pressure Pumps Cutaways models are carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, fuel system, cross-sections, lubricating circuits, cooling system etc.

Technical Specification :

Operated: Manually
Dimensions: 25 x 20 x 20(h) cm
Net Weight: 5 kg
Gross Weight: 8 kg
Equipment Layout: On Base.



Laboratory Instrument India