



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Power Assisted Steering With Electrically Controlled  
Hydraulic Pump Cutaway

**Product Code :**  
LBNY-0005-121000365



**Description :**

The Power Assisted steering consists of an electro-hydraulic system allowing power saving because the pump is driven by an electric motor.

**Technical Specification :**

Operated: Manually (Through a Hand Wheel)  
Dimensions: 42 x 64 x 40(h) cm  
Net Weight: 8 kg  
Gross Weight: 14 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**