



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Assisted Steering With Electrically Controlled
Hydraulic Pump Cutaway

Product Code :
LBNY-0005-121000365



Description :

The Power Assisted steering consists of an electro-hydraulic system allowing power saving because the pump is driven by an electric motor.

Technical Specification :

Operated: Manually (Through a Hand Wheel)
Dimensions: 42 x 64 x 40(h) cm
Net Weight: 8 kg
Gross Weight: 14 kg
Equipment Layout: On Base.



Laboratory Instrument India