



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rotary Injection Pump Cutaway

Product Code :
LBNY-0005-121000368



Description :

The distributor during the gas oil pumping and distribution phase, the transfer pump, the speed governor, the automatic advance regulator and the fuel circuit are clearly displayed. The most popular rotating pump Bosch VE has been sectioned for training purposes to show all its operating parts. It is supplied complete with an indirect injector.

Technical Specification :

Operated: Manually
Dimensions: 30 x 15 x 30(h) cm
Net Weight: 5 kg
Gross Weight: 8 kg
Equipment Layout: On Base.



Laboratory Instrument India