



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Multi process MMA, TAGS, MAGS

**Product Code :**  
LBNY-0004-740005



**Description :**

Workshop Multi process MMA, TAGS, MAGS

**Technical Specification :**

Multi process MMA, TAGS, MAGS

Output Current: DC, 250 -350 Amp

Welds mild steel, weld steel, stainless steel, and cast iron with a powerful 250Amp IGBT double voltage. Unique control mode improves welding machine reliability.

Portable Design: Features a smart digital display and alarm LED to read output welding current. Equipped with adjustable shoulder strap for convenient carry. Built-in anti-stick, dig-arc force, hot start, and other functions make it easy to operate and safe for beginners.

Digital Control Safety: Select welding electrode diameter and the recommended welding current will be automatically switched. Adjust welding current according to your welding habit for added safety.

Workshop Multi process MMA, TAGS, MAGS, Workshop Multi process MMA, TAGS, MAGS Equipments, Workshop Multi process MMA, TAGS, MAGS Tools, Workshop Multi process MMA, TAGS, MAGS Tool Kits, Workshop Multi process MMA, TAGS, MAGS Manufacturers, Workshop Multi process MMA, TAGS, MAGS Suppliers from India, China, Kenya

---

# Laboratory Instrument India

