



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Robotised Tiptronic Gearbox Cutaway

**Product Code :**  
LBNY-0005-121000371



**Description :**

Electronically-assisted manual gearbox in which the drivers command is translated into motion of the clutch sleeve timers by means of levers similar to those of the mechanical gearbox, operated by electromagnetic, electro-hydraulic or electro-pneumatic actuators.

**Technical Specification :**

Operated: Manually (Through a Crank Handle)  
Dimensions: 60 x 50 x 80(h) cm  
Net Weight: 35 kg  
Gross Weight: 55 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**