



Email : sales@laboratoryinstrumentindia.com

Product Name :

HGV Automatic Five Speed Gearbox Cutaway

Product Code :

LBNY-0005-121000372



www.laboratoryinstrumentindia.com

Description :

Sectioned heavy truck gearbox.
Composed by:
Hydrodynamic retarder
Rotating multi-disc clutches
Torque converter with lock-up clutch
Electro-valves
Fixed-position multi-disc brakes
Oil cooler with oil-water exchanger.

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 100 x 70 x 85(h) cm
Net Weight: 295 kg
Gross Weight: 350 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India

