



Email : sales@laboratoryinstrumentindia.com

Product Name :
RWD Three Speed Automatic Transmission

Product Code :
LBNY-0005-121000373



Description :

3 Forward speeds and reverse
Torque converter
Casing
Oil pump
Stationary plate clutch
Hydraulic circuit valve
Centrifugal regulator
Rotary plate clutch
Planetary gear train.

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 80 x 60 x 80(h) cm
Net Weight: 50 kg
Gross Weight: 70 kg
Equipment Layout: On Stand with Wheels.

Laboratory Instrument India

