




Email : sales@laboratoryinstrumentindia.com

Product Name :
Electric Power Assisted Steering With Suspension Cutaway

Product Code :
LBNY-0005-121000376



Description :

An instrument cluster and two displays show all information concerning the operation. Electrical power steering trainer in working conditions with fully operating McPherson suspension, for the study of the steering system. Any information used by the steering (car speed, alternator, city, steering effort) can be varied by the user.

Technical Specification :

Operated: Electrically (Motor 230 V)
Dimensions: 160 x 130 x 150(h) cm
Net Weight: 150 kg
Gross Weight: 230 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India