



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Steering Box With Hourglass Screw Cutaway

**Product Code :**  
LBNY-0005-121000377



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Steering Box With Hourglass Screw Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

**Technical Specification :**

Operated: Manually (Through a Hand Wheel)  
Dimensions: 30 x 40 x 40(h) cm  
Net Weight: 7 kg  
Gross Weight: 10 kg  
Equipment Layout: On Base.



**Laboratory Instrument India**