



Email : sales@laboratoryinstrumentindia.com

Product Name :
Steering Box With Hourglass Screw Cutaway

Product Code :
LBNY-0005-121000377



Description :

Many parts have been chromium-plated and galvanized for a longer life. Steering Box With Hourglass Screw Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections.

Technical Specification :

Operated: Manually (Through a Hand Wheel)
Dimensions: 30 x 40 x 40(h) cm
Net Weight: 7 kg
Gross Weight: 10 kg
Equipment Layout: On Base.



Laboratory Instrument India