



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Gas Regulators for MAG, TAGS Plants

Product Code :
LBNY-0004-740006



Description :

Workshop Gas Regulators for MAG, TAGS Plants

Technical Specification :

Gas Regulators for MAG, TAGS Plants

Nominal Flow Rate: 20L/min

Applicable Gas: Argon gas

Measuring Range: 0-25 (MPa)

Rated Outlet Pressure: 0-25 (MPa)

Rated flow rate 0-200: (m³ / h)

Ambient temperature 0-30

Workshop Gas Regulators for MAG, TAGS Plants, Workshop Gas Regulators for MAG, TAGS Plants Equipments, Workshop Gas Regulators for MAG, TAGS Plants Tools, Workshop Gas Regulators for MAG, TAGS Plants Tool Kits, Workshop Gas Regulators for MAG, TAGS Plants Manufacturers, Workshop Gas Regulators for MAG, TAGS Plants Suppliers from India, China, Kenya



Laboratory Instrument India

