



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radial Piston Injection Pump With Electronic Adjuster
Cutaway

Product Code :
LBNY-0005-121000383



Description :

Electric fan
Distributor plunger
Vane pump
Injection phase
Electronic control unit, etc.

Technical Specification :

Operated: Manually
Dimensions: 40 x 20 x 30(h) cm
Net Weight: 9 kg
Gross Weight: 15 kg
Equipment Layout: On Base.



Laboratory Instrument India

