



Email : sales@laboratoryinstrumentindia.com

Product Name :
HGV Sixteen Speed Manual Gearbox Cutaway

Product Code :
LBNY-0005-121000388



Description :

The over-gear allows to divide each gear into slow or fast obtaining 16 forward gears which can be inserted and geared down in sequence. The gearbox is composed of a central box containing 4 forward speeds gearings and 2 reverse speeds gearings, epicyclic unit for selecting the speed-gears mounted on the base and over-gear on top.

Technical Specification :

Operated: Manually (Through a Crank Handle)
Dimensions: 90 x 110 x 130(h) cm
Net Weight: 280 kg
Gross Weight: 380 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India