




Email : sales@laboratoryinstrumentindia.com

Product Name :
CAV DPC Injection Pump Cutaway

Product Code :
LBNY-0005-121000393



Description :

It has been cross sectioned according to the same criteria to show its main parts. CAV DPC Injection Pump Cutaway an evolution of the CAV DPA pump. It is provided with an indirect injector.

Technical Specification :

Operated: Manually
Dimensions: 30 x 20 x 20(h) cm
Net Weight: 6 kg
Gross Weight: 9 kg
Equipment Layout: On Base.



Laboratory Instrument India