



Email : sales@laboratoryinstrumentindia.com

Product Name :

Continuously Variable Transmission CVT Gearbox Cutaway

Product Code :

LBNY-0005-121000395



Description :

There are two V-belt pulleys that are split perpendicular to their axes of rotation, with a V-belt running between them. Special gearbox that can change continuously through an infinite number of effective gear ratios between maximum and minimum values.

Technical Specification :

Electro-magnetic clutch
Hydraulic command unit
Roller belt
Primary pulley
Secondary pulley
Final reduction unit
Forward and reverse insertion lever.

Other Specifications:

Operated: Manually (Through a Crank Handle)
Dimensions: 60 x 50 x 80(h) cm
Net Weight: 45 kg
Gross Weight: 75 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India