



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Continuously Variable Transmission CVT Gearbox Cutaway

**Product Code :**

LBNY-0005-121000395



**Description :**

There are two V-belt pulleys that are split perpendicular to their axes of rotation, with a V-belt running between them. Special gearbox that can change continuously through an infinite number of effective gear ratios between maximum and minimum values.

**Technical Specification :**

Electro-magnetic clutch  
Hydraulic command unit  
Roller belt  
Primary pulley  
Secondary pulley  
Final reduction unit  
Forward and reverse insertion lever.

**Other Specifications:**

Operated: Manually (Through a Crank Handle)  
Dimensions: 60 x 50 x 80(h) cm  
Net Weight: 45 kg  
Gross Weight: 75 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**