



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Workshop Welding Magnets

**Product Code :**  
LBNY-0004-740008



**Description :**

Workshop Welding Magnets

**Technical Specification :**

Welding Magnets  
30° / 60° edge, with the standard 45° and right angles  
Large centre holes for easier handling  
Weight: 0.5kg  
Two riveted holes for linking magnets together  
Spec: 45° x 60° x 90°  
Dimensions: 140mm x 112mm x 19mm  
Magnet Force: 25kg x 27kg x 37kg

Workshop Welding Magnets, Workshop Welding Magnets Equipments, Workshop Welding Magnets Tools, Workshop Welding Magnets Tool Kits, Workshop Welding Magnets Manufacturers, Workshop Welding Magnets Suppliers from India, China, Kenya



**Laboratory Instrument India**

