



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Welding Magnets

Product Code :
LBNY-0004-740008



Description :

Workshop Welding Magnets

Technical Specification :

Welding Magnets
30° / 60° edge, with the standard 45° and right angles
Large centre holes for easier handling
Weight: 0.5kg
Two riveted holes for linking magnets together
Spec: 45° x 60° x 90°
Dimensions: 140mm x 112mm x 19mm
Magnet Force: 25kg x 27kg x 37kg

Workshop Welding Magnets, Workshop Welding Magnets Equipments, Workshop Welding Magnets Tools, Workshop Welding Magnets Tool Kits, Workshop Welding Magnets Manufacturers, Workshop Welding Magnets Suppliers from India, China, Kenya



Laboratory Instrument India

