



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pick Up Baler Model Cutaway

**Product Code :**  
LBNY-0005-121000401



www.laboratoryinstrumentindia.com

**Description :**

The model clearly shows its main elements. Moreover, it has been sectioned in order to show the working principles of the internal parts. This is a model of a Pick-Up Baler which is a dragged agricultural machine. Pick-Up Baler Model Cutaway is used for the collection of the forage which is also pressed into balers.

Pick-Up Baler Model Cutaway is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts.

Pick-Up Baler Model Cutaway is composed by 4 main parts:

The woven-wire

The collecting device

A tying system

A piston system

**Technical Specification :**

Operated: Electrically (Motor 230 V)

Dimensions: 45 x 80 x 30(h) cm

Net Weight: 15 kg

Gross Weight: 25 kg

Equipment Layout: On Base.

---

# Laboratory Instrument India

