



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radial Engine Cutaway

Product Code :
LBNY-0005-121000404



Description :

As it is mechanically simple and sturdy, it is used for tanks, hovercrafts, etc. This kind of engine was largely used in aeronautic before the introduction of reaction engines.

Technical Specification :

Driving shaft with integral master rod and moving connection rods
Double ignition
Air cooling.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 150 x 140 x 180(h) cm
Net Weight: 700 kg
Gross Weight: 850 kg
Equipment Layout: On Stand with Wheels.



Laboratory Instrument India

