



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Radial Engine Cutaway

**Product Code :**  
LBNY-0005-121000404



**Description :**

As it is mechanically simple and sturdy, it is used for tanks, hovercrafts, etc. This kind of engine was largely used in aeronautic before the introduction of reaction engines.

**Technical Specification :**

Driving shaft with integral master rod and moving connection rods  
Double ignition  
Air cooling.

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 150 x 140 x 180(h) cm  
Net Weight: 700 kg  
Gross Weight: 850 kg  
Equipment Layout: On Stand with Wheels.



**Laboratory Instrument India**

