



Email : sales@laboratoryinstrumentindia.com

Product Name :
Atomiser Cutaway

Product Code :
LBNY-0005-121000405



Description :

Many parts have been chromium-plated and galvanized for a longer life. Atomiser Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits.

Technical Specification :

Hydraulic stirrer
Adjustable nozzles
Cardan shaft
Radial piston pump
Fibreglass tank with filter
Overdrive with fan
Control unit.

Other Specifications:

Operated: Electrically (Motor 230 V)
Dimensions: 220 x 110 x 130(h) cm
Net Weight: 220 kg
Gross Weight: 330 kg
Equipment Layout: Standalone.



Laboratory Instrument India