



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Atomiser Cutaway

**Product Code :**  
LBNY-0005-121000405



**Description :**

Many parts have been chromium-plated and galvanized for a longer life. Atomiser Cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits.

**Technical Specification :**

Hydraulic stirrer  
Adjustable nozzles  
Cardan shaft  
Radial piston pump  
Fibreglass tank with filter  
Overdrive with fan  
Control unit.

**Other Specifications:**

Operated: Electrically (Motor 230 V)  
Dimensions: 220 x 110 x 130(h) cm  
Net Weight: 220 kg  
Gross Weight: 330 kg  
Equipment Layout: Standalone.



**Laboratory Instrument India**